

### Attempt 1

1. Manual Delineation [Measure mCNV area by measuring pixels in contour]
  - i) Selection using Polygon Selection Tool
  - ii) Edit->Clear Outside
  - iii) Image->Type->8-bit
2. Smooth Image (Gaussian Kernel)
  - i) Process->Filter->Gaussian Blur (Radius =1)
3. Binary image (Frangi vesselness filter)
  - i) Process->Filters->Frangi Vesselness (Apply Gaussian to all scales = True, Spacing = 1,1, Scale = 3,5)
  - ii) Image->Type->8 bit
4. Binary Image (local adaptiveness thresholding) [Measure vessel area, fractal dimension, vessel lengths]
  - For "Before"
    - i) Image->Adjust->Auto Threshold (Method = Mean)
    - ii) Process->Binary->Close
  - For "After"
    - i. Image->Adjust->Auto Local Threshold (Method = Mean, Radius = 15)
5. Tagged Skeleton Image [Calculate number of vessels]
  - i) Process->Binary->Skeletonise
  - ii) LUT->Invert LUT

### Attempt 2